

## Product datasheet for **MP206842**

### Iqce Mouse qPCR Primer Pair (NM\_028833)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74239
Forward Sequence:	AGAGCATCTCCAAGGAACCGTG
Reverse Sequence:	TCCTTCTCCAGGTCGATCTTGG
ACCN:	<u><a href="#">BC059223</a></u> , <u><a href="#">NM_028833</a></u> , <u><a href="#">NM_028833.1</a></u> , <u><a href="#">NM_028833.2</a></u> , <u><a href="#">NM_028833.3</a></u> , <u><a href="#">BC042766</a></u> , <u><a href="#">BC065128</a></u>
UniProt ID:	<u><a href="#">Q6PCQ0</a></u>
Synonyms:	1700028P05Rik; mKIAA1023
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)