

Product datasheet for **MP206899**

Kdm4b Mouse qPCR Primer Pair (NM_172132)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	193796
Forward Sequence:	TGGTGTCCATCACCTGCCTCAA
Reverse Sequence:	CATTGCGGTTCTTGGTGATGACG
ACCN:	<u>BC007145, NM_172132, NM_172132.1, NM_172132.2, BB213288, BC005480, NM_172132.3</u>
UniProt ID:	<u>Q91VY5</u>
Synonyms:	4732474L06Rik; Jmjd2b; mKIAA0876
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)