

Product datasheet for **MP206923**

Jph4 Mouse qPCR Primer Pair (NM_177049)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319984
Forward Sequence:	CAAGCTGATAGCCCAGGATCTAC
Reverse Sequence:	GCTGTCTCATCCAAAGGCTCT
ACCN:	BC052049 , BC052376 , BC052722 , NM_177049 , NM_177049.1 , NM_177049.2 , NM_177049.3 , NM_177049.4 , NM_177049.5 , BC052049.1
UniProt ID:	Q80WT0
Synonyms:	9330157P13Rik; JP-4; JPHL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)