

Product datasheet for **MP206955**

Kpna2 Mouse qPCR Primer Pair (NM_010655)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16647
Forward Sequence:	CAAACACGGTGCGATTGACCCA
Reverse Sequence:	GAAGAATCTGCTCAACGGCGTC
ACCN:	<u>BC003274</u> , <u>BC006720</u> , <u>BC082280</u> , <u>BC094011</u> , <u>NM_010655</u> , <u>NM_010655.1</u> , <u>NM_010655.2</u> , <u>NM_010655.3</u> , <u>BC003274.1</u> , <u>BC086789</u> , <u>BC106110</u>
UniProt ID:	<u>P52293</u>
Synonyms:	2410044B12Rik; IPOA1; PTAC58; Rch1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)