

Product datasheet for **MP206964**

Kazald1 Mouse qPCR Primer Pair (NM_178929)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107250
Forward Sequence:	TGGACGCTAACCTCACTGTGGT
Reverse Sequence:	GGCAAACACCTCACAGCCAAAG
ACCN:	BC115642 , NM_178929 , NM_178929.1 , NM_178929.2 , NM_178929.3 , NM_178929.4 , BB618180 , BB871162 , BC115643 , BX636599 , NM_178929.5
UniProt ID:	Q8BJ66
Synonyms:	AI842353; Bono1; IGFBP-rP10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)