

Product datasheet for **MP206974**

Klhl42 Mouse qPCR Primer Pair (NM_001081237)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	232539
Forward Sequence:	AGTTCTCTGCCTGCGAGTGCAA
Reverse Sequence:	GTCACAGGTCTCAAACAGCGTC
ACCN:	<u>NM_001081237</u> , <u>NM_001081237.1</u> , <u>BC151076</u> , <u>BC011326</u> , <u>BC030345</u> , <u>BC037084</u> , <u>BC147633</u> , <u>BC147636</u> , <u>BC151064</u> , <u>BC151065</u> , <u>BC151068</u>
UniProt ID:	<u>Q8BFQ9</u>
Synonyms:	A330079N23Rik; AI450355; C230080I20Rik; Ctb9; Klhdc5; mKIAA1340
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)