

## Product datasheet for **MP207056**

### Krt77 Mouse qPCR Primer Pair (NM\_001003667)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Gene ID:          | 406220                                                                                                                                                                                                                                                                                                                                            |
| Forward Sequence: | GCTCAACCGTAACATCCAGAGG                                                                                                                                                                                                                                                                                                                            |
| Reverse Sequence: | TTGGCGTTCTGGAGAGCCCTTT                                                                                                                                                                                                                                                                                                                            |
| ACCN:             | <u><a href="#">BC107395</a></u> , <u><a href="#">BC107396</a></u> , <u><a href="#">NM_001003667</a></u> , <u><a href="#">NM_001003667.1</a></u> , <u><a href="#">BC141259</a></u> , <u><a href="#">BK003993</a></u> , <u><a href="#">BK005256</a></u> ,<br><u><a href="#">NM_001003667.2</a></u>                                                  |
| UniProt ID:       | <u><a href="#">Q6IFZ6</a></u>                                                                                                                                                                                                                                                                                                                     |
| Synonyms:         | 4732484G22Rik; K1emb; Krt239; mK1b                                                                                                                                                                                                                                                                                                                |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                       |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | Store at -20 °C.                                                                                                                                                                                                                                                                                                                                  |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                      |

[View online »](#)