

Product datasheet for **MP207100**

Khdrbs2 Mouse qPCR Primer Pair (NM_133235)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	170771
Forward Sequence:	GCTTGGACCAAGAGGAACTCC
Reverse Sequence:	TCACTCAGGTGGGCATACTGG
ACCN:	<u>BC132117, BC132119, NM_133235, NM_133235.1, NM_133235.2, NM_133235.3, BC066814</u>
UniProt ID:	<u>Q9WU01</u>
Synonyms:	6330586C16Rik; mSLM-1; Slim1; SLM; SIm-1; SIm1; Tg(LRRK2*R1441G)135Cjli; TG-RP135
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)