

Product datasheet for **MP207133**

Klra1 Mouse qPCR Primer Pair (NM_016659)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16627
Forward Sequence:	AGACCTGAGGAGACTAAAGGGC
Reverse Sequence:	GCCAACACTGAAACAGCTACCAG
ACCN:	<u>NM_016659</u> , <u>NM_016659.1</u> , <u>NM_016659.2</u> , <u>NM_016659.3</u> , <u>BC140968</u> , <u>BC144929</u>
Synonyms:	A1; CH29-493D4.3; Klra22; Ly49a; Ly49o<129>; Ly49v
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)