

Product datasheet for **MP207140**

Klra3 Mouse qPCR Primer Pair (NM_010648)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16634
Forward Sequence:	GAGAACAGGACAGATGGGACAG
Reverse Sequence:	GCAGTTCGCTTTACATCCACTCC
ACCN:	<u>BC099867</u> , <u>BC099878</u> , <u>BC099879</u> , <u>BC099880</u> , <u>BC101952</u> , <u>BC103543</u> , <u>BC103544</u> , <u>BC103545</u> , <u>BC132524</u> , <u>BC132526</u> , <u>NM_010648</u> , <u>NM_010648.1</u> , <u>NM_010648.2</u> , <u>NM_010648.3</u> , <u>BC141217</u> , <u>BC145192</u> , <u>BC145193</u> , <u>BX515458</u>
UniProt ID:	<u>Q64329</u>
Synonyms:	5E6; Ly49c; Nk-2; NK-2.1; Nk2; Nk2.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)