

## Product datasheet for **MP207150**

### Kdr Mouse qPCR Primer Pair (NM\_010612)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16542
Forward Sequence:	CGAGACCATTGAAGTGACTTGCC
Reverse Sequence:	TTCCTCACCTGCGGATAGTCA
ACCN:	<a href="#">BC020530</a> , <a href="#">NM_010612</a> , <a href="#">NM_010612.1</a> , <a href="#">NM_010612.2</a> , <a href="#">NM_010612.3</a>
UniProt ID:	<a href="#">P35918</a>
Synonyms:	6130401C07; Flk-1; Flk1; Krd-1; Ly73; sVEGFR-2; VEGFR-2; VEGFR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)