

Product datasheet for **MP207205**

Kank1 Mouse qPCR Primer Pair (NM_181404)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107351
Forward Sequence:	CAGTTCATTGGCATTAAACGGCGG
Reverse Sequence:	AGTGTCGTGGTCCTTCTCTTCC
ACCN:	<u>BC079563</u> , <u>NM_181404</u> , <u>NM_181404.1</u> , <u>NM_181404.2</u> , <u>NM_181404.3</u> , <u>NM_181404.4</u> , <u>NM_181404.5</u> , <u>NM_181404.6</u> , <u>BC037712</u> , <u>BC050250</u> , <u>BC058673</u> , <u>BC152400</u>
Synonyms:	A930031B09Rik; Ankrd15; AU015049; AW121052; D330024H06Rik; mKIAA0172
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)