

Product datasheet for **MP207231**

Poglut1 Mouse qPCR Primer Pair (NM_172380)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	224143
Forward Sequence:	GAAAGCCAGCCGCTAAGGATGT
Reverse Sequence:	TGGAAGACCAGTGAACCGCACA
ACCN:	<u>BC026809, NM_172380, NM_172380.1, NM_172380.2, NM_172380.3, NM_172380.4</u>
UniProt ID:	<u>Q8BYB9</u>
Synonyms:	9630046K23Rik; Clp46; Ktelc; Ktelc1; Ru; Rumi; w; wsnp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)