

Product datasheet for **MP207249**

Larp4b Mouse qPCR Primer Pair (NM_172585)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	217980
Forward Sequence:	GCCAAGTAGTCTTGAATCCTCTC
Reverse Sequence:	TGGTGTGGAGACTGCGTCTGA
ACCN:	<u>BC078640</u> , <u>NM_172585</u> , <u>NM_172585.1</u> , <u>NM_172585.2</u> , <u>NM_172585.3</u> , <u>BC051459</u> , <u>BC059218</u> , <u>BC086626</u> , <u>BU056450</u>
UniProt ID:	<u>Q6A0A2</u>
Synonyms:	A630096F19; AI256361; D13Wsu64e; Larp5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)