

Product datasheet for **MP207278**

Lamb1 Mouse qPCR Primer Pair (NM_008482)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16777
Forward Sequence:	GAACTACACGGTGAGGTTGGAG
Reverse Sequence:	GCCAACAGTGAAGATGTCCAGC
ACCN:	<u>NM_008482</u> , <u>NM_008482.1</u> , <u>NM_008482.2</u> , <u>BC150809</u> , <u>BC032276</u> , <u>BC121790</u>
Synonyms:	C77966; C80098; C81607; D130003D08Rik; Lamb-1; Lamb1-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)