

Product datasheet for **MP207282**

Lama4 Mouse qPCR Primer Pair (NM_010681)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16775
Forward Sequence:	CAGTTTGTCTCTACCTCGGAAG
Reverse Sequence:	CTCACAGGCTTGGATCCAGGA
ACCN:	BC115942 , NM_010681 , NM_010681.1 , NM_010681.2 , NM_010681.3 , NM_010681.4 , BC008533 , BC038526 , BC056176
UniProt ID:	P97927
Synonyms:	MGC141567
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)