

## Product datasheet for **MP207317**

### Adgrl1 Mouse qPCR Primer Pair (NM\_181039)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	330814
Forward Sequence:	CCATCAAGCAGAACAGCCGCAA
Reverse Sequence:	GCTTCACTGTGGCATTCTCCGT
ACCN:	<u><a href="#">BC085138</a></u> , <u><a href="#">NM_181039</a></u> , <u><a href="#">NM_181039.1</a></u> , <u><a href="#">NM_181039.2</a></u> , <u><a href="#">BC026388</a></u> , <u><a href="#">BC055793</a></u> , <u><a href="#">BC078414</a></u>
UniProt ID:	<u><a href="#">Q80TR1</a></u>
Synonyms:	2900070I05Rik; AI182398; CLIBA; Lec2; Lphn1; mKIAA0821
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)