

Product datasheet for **MP207328**

Lgals4 Mouse qPCR Primer Pair (NM_010706)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16855
Forward Sequence:	CTTCAGTCCATCAACTTCCTCGG
Reverse Sequence:	GGCAAGGTGTTTCATCTGTGGAG
ACCN:	BC011236 , BC021632 , BC030297 , BC094008 , NM_010706 , NM_010706.1 , NM_010706.2 , BC141106 , BC145387
UniProt ID:	Q8K419
Synonyms:	gal-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)