

Product datasheet for **MP207331**

Lgals8 Mouse qPCR Primer Pair (NM_018886)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56048
Forward Sequence:	GAGGAGATCACCTACGACATGC
Reverse Sequence:	CGTACAGCAGAACATGCCTTCC
ACCN:	<u>BC040243</u> , <u>NM_018886</u> , <u>NM_018886.1</u> , <u>NM_018886.2</u> , <u>NM_018886.3</u> , <u>NM_018886.4</u> , <u>NM_018886.5</u> , <u>BQ033229</u>
UniProt ID:	<u>Q9JL15</u>
Synonyms:	1200015E08Rik; AI326142; D13ErtD524e; Lgals-8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)