

Product datasheet for **MP207334**

Lgals3bp Mouse qPCR Primer Pair (NM_011150)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19039
Forward Sequence:	ATGCCTGTCGTCAGAGACTTCC
Reverse Sequence:	GGCAAAAAGCCGTCCACAGTAG
Accession No:	<u>BC090658</u> , <u>NM_011150</u>
Synonyms:	90K; CyCAP; MAC-2BP; Ppicap; Tango10b
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)