

Product datasheet for **MP207339**

Lefty1 Mouse qPCR Primer Pair (NM_010094)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	13590
Forward Sequence:	AGTCCTGGACAAGGCTGATGTG
Reverse Sequence:	GAGGTCTCTGACACCAGGAACC
ACCN:	<u>BC050221, NM_010094, NM_010094.1, NM_010094.2, NM_010094.3, NM_010094.4</u>
UniProt ID:	<u>Q64280</u>
Synonyms:	AI450052; Eba; Ebaf; Left; Leftb; lefty; lefty-; lefty-1; St; Stra3; Tgfb4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)