

Product datasheet for **MP207349**

Lmtk2 Mouse qPCR Primer Pair (NM_001081109)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	231876
Forward Sequence:	AGTCATTCGGGAGAGGGACATG
Reverse Sequence:	CCTCAGGTAAGTTAGCAGACGG
ACCN:	<u>NM_001081109</u> , <u>NM_001081109.1</u> , <u>BC058546</u> , <u>BC058653</u> , <u>NM_001081109.2</u>
UniProt ID:	<u>Q3TYD6</u>
Synonyms:	2900041G10Rik; A330101P12Rik; aatyk2; BREK; cprk; KPI-2; KPI2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)