

Product datasheet for **MP207385**

Carmil1 Mouse qPCR Primer Pair (NM_026825)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68732
Forward Sequence:	GCTGTCATCCAGTCCCATTTCAG
Reverse Sequence:	TCCTGTGCGATGTGGGCTTTTC
ACCN:	<u>NM_026825</u> , <u>NM_026825.1</u> , <u>NM_026825.2</u> , <u>NM_026825.3</u> , <u>BC012229</u> , <u>BC150819</u> , <u>BC167257</u>
UniProt ID:	<u>Q6EDY6</u>
Synonyms:	1110037D04Rik; AI425970; CARMIL; CARML1; D130057M20; Lrrc16; Lrrc16a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)