

Product datasheet for **MP207387**

Lrrc18 Mouse qPCR Primer Pair (NM_026253)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67580
Forward Sequence:	GGAGAAAACGCCTTGACCTGAG
Reverse Sequence:	CAGGTTCTGGAACCTGGCGATG
ACCN:	<u>BC061072</u> , <u>NM_026253</u> , <u>NM_026253.1</u> , <u>NM_026253.2</u> , <u>NM_026253.3</u> , <u>NM_026253.4</u> , <u>BB017856</u> , <u>BC048444</u>
UniProt ID:	<u>Q9CQ07</u>
Synonyms:	4930442L21Rik; MTLR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)