

Product datasheet for **MP207396**

Lrrc28 Mouse qPCR Primer Pair (NM_175124)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67867
Forward Sequence:	AGGAACTCCCTCACAACCTTGC
Reverse Sequence:	CAGACACTGGAGTTTCACGAGG
ACCN:	BC086694 , NM_175124 , NM_175124.1 , NM_175124.2 , NM_175124.3 , NM_175124.4 , BC027058 , BC032971 , BC052051
UniProt ID:	Q3TX51
Synonyms:	1300004K21Rik; 2210012C09Rik; 2310058O11Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)