

Product datasheet for **MP207398**

Lrrc3 Mouse qPCR Primer Pair (NM_145152)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237387
Forward Sequence:	TCTGCTTACGGTTAGGTGCCTC
Reverse Sequence:	GGCGTCTAGTTTCAGAAGCACC
ACCN:	BC115652 , NM_145152 , NM_145152.1 , NM_145152.2 , NM_145152.3 , NM_145152.4 , BC058363 , BQ176232 , BY320829
UniProt ID:	P59034
Synonyms:	1300011L04Rik; C21orf102
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)