

Product datasheet for **MP207407**

Lrrc41 Mouse qPCR Primer Pair (NM_153521)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	230654
Forward Sequence:	GCTCTGACACTGTGCCATCTTC
Reverse Sequence:	CTGTCCAGGAGGCGGAAGATAT
ACCN:	BC022976 , NM_153521 , NM_153521.1 , NM_153521.2 , BB619031 , BC005765 , BC007169 , BC024940 , BC071258 , BQ032399 , BY033356
UniProt ID:	Q8K1C9
Synonyms:	AA409966; AW555107; D630045E04Rik; D730026A16Rik; MUF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)