

Product datasheet for **MP207413**

Lrrc47 Mouse qPCR Primer Pair (NM_201226)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	72946
Forward Sequence:	CTGGATGGGAAGGAGAACTACC
Reverse Sequence:	TGCAGGCTAGTGGCACTCGTTA
ACCN:	<u>BC094573</u> , <u>NM_201226</u> , <u>NM_201226.1</u> , <u>BC048911</u> , <u>BC138462</u> , <u>BC138463</u> , <u>BU148242</u> , <u>BU514662</u>
UniProt ID:	<u>Q505F5</u>
Synonyms:	2900010D03Rik; AA960559; mKIAA1185
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)