

Product datasheet for **MP207441**

Lrrtm4 Mouse qPCR Primer Pair (NM_178731)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	243499
Forward Sequence:	TTCCTCGTCTCTCAACCTGCG
Reverse Sequence:	CAGGCTCGATTGCTTGGATGTC
ACCN:	<u>BC037216</u> , <u>NM_178731</u> , <u>NM_178731.1</u> , <u>NM_178731.2</u> , <u>NM_178731.3</u> , <u>NM_178731.4</u> , <u>NM_178731.5</u> , <u>BC037216.1</u> , <u>BY130176</u>
UniProt ID:	<u>Q80XG9</u>
Synonyms:	7530419J18Rik; A230052N11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)