

Product datasheet for **MP207453**

Pidd1 Mouse qPCR Primer Pair (NM_022654)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57913
Forward Sequence:	ACAGAAGAGCCTCGGCAAGTCT
Reverse Sequence:	GAAAGGCACAGCAGAGGGCTTA
ACCN:	<u>NM_022654</u> , <u>NM_022654.1</u> , <u>BC145857</u> , <u>BC032202</u> , <u>BC056987</u> , <u>BC089373</u> , <u>BC099686</u> , <u>BC145314</u> , <u>BC145315</u> , <u>NM_022654.2</u>
UniProt ID:	<u>Q9ERV7</u>
Synonyms:	1200011D09Rik; AU042446; Lrdd; Pidd
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)