

Product datasheet for **MP207463**

Lgr4 Mouse qPCR Primer Pair (NM_172671)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107515
Forward Sequence:	CAAATGCCACCAGCACTGCTGA
Reverse Sequence:	CCCACTGTAAGGCGAATCATCC
ACCN:	<u>NM_172671</u> , <u>NM_172671.1</u> , <u>NM_172671.2</u> , <u>BC056637</u> , <u>BC160293</u> , <u>BE861913</u> , <u>BE987673</u> , <u>BU520309</u>
UniProt ID:	<u>A2ARI4</u>
Synonyms:	9130225G07; A330106J01Rik; A930009A08Rik; Gpr48
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)