

Product datasheet for **MP207468**

Lrch4 Mouse qPCR Primer Pair (NM_146164)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	231798
Forward Sequence:	GATGGAGACCTGGAGCAGATTG
Reverse Sequence:	GCTGCTAGATGGCTTCTCCACA
ACCN:	<u>BC012525</u> , <u>NM_146164</u> , <u>NM_146164.1</u> , <u>NM_146164.2</u> , <u>BQ924416</u> , <u>BY204138</u>
UniProt ID:	<u>Q921G6</u>
Synonyms:	2810008P14Rik; 2900069C24Rik; AI558103; LRN; LRRN4; mFLJ00248; SAP25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)