

Product datasheet for **MP207475**

Lztr1 Mouse qPCR Primer Pair (NM_025808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66863
Forward Sequence:	CGGCTTCTACAATAACCGACTGC
Reverse Sequence:	ATGTCCAGTGCCTGCGTCTTGT
ACCN:	<u>BC034400</u> , <u>NM_025808</u> , <u>NM_025808.1</u> , <u>NM_025808.2</u> , <u>NM_025808.3</u> , <u>NM_025808.4</u> , <u>BC022743</u> , <u>BF457766</u> , <u>BY005614</u> , <u>BY146849</u>
UniProt ID:	<u>Q9CQ33</u>
Synonyms:	I200003E2IRik; AI591627; AW550890; TCFL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

