

Product datasheet for **MP207481**

Cnmd Mouse qPCR Primer Pair (NM_010701)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16840
Forward Sequence:	TCCTTGAACCTCTGTGGCGACCT
Reverse Sequence:	GGAGCACTGTTTCTCAGACTTC
ACCN:	BC045152 , NM_010701 , NM_010701.1 , NM_010701.2 , NM_010701.3 , BB849764 , BB849846 , BU705377 , BY776014 , BY785008
UniProt ID:	Q9Z1F6
Synonyms:	Bricd3; ChM-I; Chmd; Lect1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)