

Product datasheet for **MP207506**

Large1 Mouse qPCR Primer Pair (NM_010687)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16795
Forward Sequence:	CGAGTTTGAGCCTTCTGCTGAC
Reverse Sequence:	GGCTTCTGCATCTGACAGGTAG
ACCN:	BC100399 , BC125487 , BC125489 , NM_010687 , NM_010687.1 , NM_010687.2 , BC061506 , BC079596 , BQ442282
UniProt ID:	Q9Z1M7
Synonyms:	BPFD#36; enr; fg; froggy; Gyltl1a; Large; Mbp-1; Mbp1; myd
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)