

Product datasheet for **MP207511**

Lims1 Mouse qPCR Primer Pair (NM_026148)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110829
Forward Sequence:	GGGTTTGTCAAGAATGCTGGCAG
Reverse Sequence:	GCACAGTTGAAGTGGTCTGGATG
ACCN:	<u>BC005621</u> , <u>NM_026148</u> , <u>NM_026148.1</u> , <u>NM_026148.2</u> , <u>NM_026148.3</u> , <u>BC067393</u> , <u>BQ554039</u> , <u>NM_026148.4</u>
UniProt ID:	<u>Q99JW4</u>
Synonyms:	2310016J22Rik; 4921524A02Rik; AI507642; AU021743; AW551584; C430041B13Rik; Lims1l; PINCH; PINCH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)