

Product datasheet for **MP207522**

Lmo4 Mouse qPCR Primer Pair (NM_010723)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16911
Forward Sequence:	ATAGCCGCTGCCTCAAGTGCTC
Reverse Sequence:	AGGCACTGCAAGCACCGCTATT
ACCN:	BC003488 , BC010278 , BC100382 , NM_010723 , NM_010723.1 , NM_010723.2 , NM_010723.3 , BC004661 , BY267947
UniProt ID:	P61969
Synonyms:	A730077C12Rik; Crp3; Etohi4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)