

Product datasheet for **MP207535**

Lnpg Mouse qPCR Primer Pair (NM_027133)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	69605
Forward Sequence:	CGCCAGAAAGAACTGTTGCTCC
Reverse Sequence:	TATCCAGAGCACCTCGTTCTCG
ACCN:	<u>BC057961</u> , <u>NM_027133</u> , <u>NM_027133.1</u> , <u>NM_027133.2</u> , <u>NM_027133.3</u> , <u>BC060153</u> , <u>BY161702</u> , <u>BY770278</u>
UniProt ID:	<u>Q7TQ95</u>
Synonyms:	2310011O18Rik; 4921514L11Rik; 9530051D01Rik; A1666268; Lnp; Lnpg1; Ul; ulnaless
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)