

Product datasheet for **MP207541**

Limk1 Mouse qPCR Primer Pair (NM_010717)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16885
Forward Sequence:	GAGAAACAGGCGAGGTGATGGT
Reverse Sequence:	CACTCCGATGAACTTGAGCACG
ACCN:	<u>NM_010717</u> , <u>NM_010717.1</u> , <u>NM_010717.2</u> , <u>NM_010717.3</u> , <u>BC063768</u> , <u>BC156778</u> , <u>BQ828670</u> , <u>BX527775</u> , <u>BY230595</u>
UniProt ID:	<u>P53668</u>
Synonyms:	KIZ-1; Limk; LIMK-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)