

Product datasheet for **MP207547**

Lin7a Mouse qPCR Primer Pair (NM_001039354)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	108030
Forward Sequence:	GGTGGCTGAAAGACATGGAGGC
Reverse Sequence:	GCTTCACACTGTCCTTAGCAGC
ACCN:	<u>BC105965</u> , <u>NM_001039354</u> , <u>NM_001039354.1</u> , <u>BC029721</u> , <u>BG288381</u> , <u>BU923246</u> , <u>BY123635</u>
UniProt ID:	<u>Q8JZS0</u>
Synonyms:	AI848705; LIN-7A; MALS-1; TIP-33; Veli; Veli1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)