

Product datasheet for **MP207549**

Lin7c Mouse qPCR Primer Pair (NM_011699)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22343
Forward Sequence:	GCCACAGAACTTCAGGCTTTGC
Reverse Sequence:	CTGCTGATGTCCACCGTCTCAT
ACCN:	<u>BC046966</u> , <u>NM_011699</u> , <u>NM_011699.1</u> , <u>NM_011699.2</u> , <u>NM_011699.3</u> , <u>BC004685</u> , <u>BC052471</u>
UniProt ID:	<u>O88952</u>
Synonyms:	9130007B12Rik; AI303698; AU019331; AW125731; D2Ert520e; LIN-7-C; LIN-7C; MALS-3; Veli3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)