

Product datasheet for **MP207556**

Lmf2 Mouse qPCR Primer Pair (NM_178919)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	105847
Forward Sequence:	CGACTTACTGCCTTCTACGCTC
Reverse Sequence:	TCAGCAGACAGGTGACGGTCAT
ACCN:	<u>BC079532</u> , <u>NM_178919</u> , <u>NM_178919.2</u> , <u>NM_178919.3</u> , <u>NM_178919.4</u> , <u>BC026448</u> , <u>BC038039</u> , <u>BC058973</u>
UniProt ID:	<u>Q8C3X8</u>
Synonyms:	AI451006; Tmem112b; Tmem153
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)