

Product datasheet for **MP207596**

Lonrf3 Mouse qPCR Primer Pair (NM_028894)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74365
Forward Sequence:	GCATGTGCCTTGGAGATCCTGT
Reverse Sequence:	CCTCTTGCCTATGCTGTCTACC
ACCN:	<u>BC120688</u> , <u>BC120690</u> , <u>NM_028894</u> , <u>BC058671</u> , <u>BY740851</u>
UniProt ID:	<u>Q9D4H7</u>
Synonyms:	4932412G04Rik; 5730439E01Rik; A830039N02Rik; AU023707; Rnf127
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)