

Product datasheet for **MP207597**

Lor Mouse qPCR Primer Pair (NM_008508)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	16939
Forward Sequence:	ACTCCTCTCAGCAGACCAGTCA
Reverse Sequence:	AGAGGAGCCACCTCCACAGCT
ACCN:	BC058223 , NM_008508 , NM_008508.1 , NM_008508.2 , BC026781
UniProt ID:	P18165
Synonyms:	AI036317; S77319
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)