

## Product datasheet for **MP207605**

### Lrp11 Mouse qPCR Primer Pair (NM\_172784)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	237253
Forward Sequence:	GACATAGCGTTGGCTTGTGATGG
Reverse Sequence:	TGGGCAGAAGTAGCCACTGTGT
ACCN:	<u>NM_172784, NM_172784.1, NM_172784.2, NM_172784.3, BC059874, NM_172784.4</u>
UniProt ID:	<u>Q8CB67</u>
Synonyms:	1700034J19Rik; 6330533B21; 9830160H19Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)