

Product datasheet for **MP207623**

Lsm2 Mouse qPCR Primer Pair (NM_030597)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27756
Forward Sequence:	GGCAAGGATGTAGTCGTGGAAC
Reverse Sequence:	AGGGTCTGTGACACTGATGTCG
ACCN:	<u>BC014288</u> , <u>BC049543</u> , <u>NM_030597</u> , <u>NM_030597.1</u> , <u>NM_030597.2</u> , <u>NM_030597.3</u> , <u>BB848823</u> , <u>BY079930</u>
UniProt ID:	<u>O35900</u>
Synonyms:	D17H6S56E-2; D17H6S56E2; Dmapl; Dmpkap; G7b; Sm-X5; SmX5; snRNP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)