

Product datasheet for **MP207633**

Luc7l2 Mouse qPCR Primer Pair (NM_138680)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	192196
Forward Sequence:	CCCGATCACATAGCAAGAATCCC
Reverse Sequence:	CGATGGCGTTTCTCAGAGACT
ACCN:	BC020058 , BC056354 , BC056383 , BC056970 , NM_138680 , NM_138680.1 , NM_138680.2 , NM_138680.3 , BC020058.1 , BC002314 , BC016229 , BC040793 , BC054449 , BQ128768
UniProt ID:	Q7TNC4
Synonyms:	4930471C18Rik; AA522013; AU015269; CGI59; Cgi74
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)