

Product datasheet for **MP207655**

Ly6g6c Mouse qPCR Primer Pair (NM_023463)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68468
Forward Sequence:	CACTCCTGCTACAAAGTCCCTG
Reverse Sequence:	CCGACAAGGCTCTTCTGGTGTG
ACCN:	<u>BC116365</u> , <u>BC116366</u> , <u>NM_023463</u> , <u>NM_023463.1</u> , <u>NM_023463.2</u> , <u>NM_023463.3</u>
UniProt ID:	<u>Q9Z1Q4</u>
Synonyms:	1110003M04Rik; AU016360; G6c; NG24
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)