

Product datasheet for **MP207656**

Ly6g6d Mouse qPCR Primer Pair (NM_033478)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114654
Forward Sequence:	CAGCAACTCCTGCAAGCAGACA
Reverse Sequence:	CTGGGTAGTGATGGCTCAAGGT
ACCN:	<u>BC131988</u> , <u>BC131990</u> , <u>NM_033478</u> , <u>NM_033478.1</u> , <u>NM_033478.2</u> , <u>BB737333</u> , <u>BM933266</u> , <u>NM_033478.3</u>
UniProt ID:	<u>Q9Z1Q3</u>
Synonyms:	A930024F17Rik; G6d; G6f; g6f-ly6g6d; MEGT1; NG25; NG32
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)