

Product datasheet for **MP207676**

Lyrm1 Mouse qPCR Primer Pair (NM_029610)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73919
Forward Sequence:	CTAAAGGAAGCCAGAACACTGTTC
Reverse Sequence:	GTGAGACCCATTGGAGGCAGAT
ACCN:	BC011339 , NM_029610 , NM_029610.1 , NM_029610.2 , NM_029610.3 , BU757765 , BY464772 , BY752674
UniProt ID:	Q9CQB7
Synonyms:	1110065L10Rik; 2310004B22Rik; 4930404J24Rik
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)